



# 2024 BROCHURE




**MS Pipette**

*Innovative Life Sciences Tools*

# MS Pipette



- Light weight
- Soft spring system for extremely low pipetting forces
- Controlled volume setting to prevent accidental volume changes
- Different color coded push button for different volume
- Contoured shape: fits either small or large hands
- Large pushbutton, rounded and freely rotating
- Finger hook: takes the weight, for a more relaxed grip
- Easy on site calibration (calibration key included)
- High accuracy and precision
- A unique serial number
- UV resistance even under prolonged exposure
- Fully autoclavable (121°C/0.1MPa/20 min)



3-position tip ejector button (left/right handed users)



## Ordering Information

Cat. No	Description
MSP-2	MS pipette, variable volume 0.2~2µl
MSP-10	MS pipette, variable volume 1~10µl
MSP-20	MS pipette, variable volume 2~20µl
MSP-100	MS pipette, variable volume 10~100µl
MSP-200	MS pipette, variable volume 20~200µl
MSP-1000	MS pipette, variable volume 100~1000µl
MSP-5000	MS pipette, variable volume 500~5000µl
MSP-10000	MS pipette, variable volume 1000~10000µl
MSP-8X10	MS 8-ch pipette, variable volume 0.5~10µl
MSP-8X20	MS 8-ch pipette, variable volume 2~20µl
MSP-8X200	MS 8-ch pipette, variable volume 20~200µl
MSP-8X300	MS 8-ch pipette, variable volume 20~300µl
MSP-12X10	MS 12-ch pipette, variable volume 0.5~10µl
MSP-12X20	MS 12-ch pipette, variable volume 2~20µl
MSP-12X200	MS 12-ch pipette, variable volume 20~200µl
MSP-12X300	MS 12-ch pipette, variable volume 20~300µl
MSP-ST03	Pipette Stand



**Taiwan Office**

No. 156, Sec. 1, Guoji Rd., Taoyuan Dist.,  
Taoyuan City 33061, Taiwan  
T/+886-3-3762878  
F/+886-3-3761310

**US Office**

19959 Sea Gull Way  
Saratoga, CA 95070  
U.S.A.  
T/ +1-408-366-9866  
F/ +1-408-446-1107

**Shanghai Office**

Room 612, International business exhibition center,  
9300 Hunan Road, Pudong, Shanghai, China  
National toll-free No.:400-823-9177  
T/ +86-21-50795277  
F/ +86-21-50795277

**India Office**

D.No.7-143, 2nd Floor,  
St.No.2,Nagendra Nagar,  
Habsiguda, Hyderabad-500007.  
India  
T/ +91-40-27001515  
T/ +91-40-27001586

